

# Low vacuum drainage set with collection bag



## DRENOMED Closet and continuous low vacuum drain set

- Suction bellow 300 ml capacity, with double non return valve
- PVC connecting tube length 110 cm, with connector for PVC and silicone drains
- Collection bag 500 ml capacity, with non return valve, antibacterial air vent
- Two non return valves and closed circuit to avoid risk of contamination for patient and operator
- Double packaged, sterile.

### SET WITH CONNECTING TUBE

Code	Description	Q.ty/box	Connector
MS/115V			universal type for PVC drains /size from Ch. 06 to 18)
MS/115VD	Closed and continuous low vacuum drain set with collection bag 500 ml and connecting tube	30	for silicone drains (flat, channelled S-Line or straight type size from Ch. 12 to 18)
MS/115VM			for silicone straight drains (size from Ch. 19 to 24)






### SET WITH CONNECTING TUBE AND PVC DRAIN / NEEDLE

Code	Description	Q.ty/box	Drain type
MS/116V/Ch	Closed and continuous low vacuum drain set with collection bag 500 ml and connecting tube and PVC drain (size from Ch. 06 to 18)		
MS/116VY/Ch	Closed and continuous low vacuum drain set with collection bag 500 ml, "Y" connecting tube and two PVC drains (size from Ch. 06 to 18)	30	
MS/117V/Ch	Closed and continuous low vacuum drain set with collection bag 500 ml, connecting tube and PVC drain with needle (size from Ch. 06 to 18)		
MS/117VY/Ch	Closed and continuous low vacuum drain set with collection bag 500 ml, "Y" connecting tube and two PVC drains with needle (size from Ch. 06 to 18)		


Note: replace "Ch" in the code with the drain size (Ch. 06, 08, 10, 12, 14, 16 or 18).

## SET WITH CONNECTING TUBE AND SILICONE DRAIN / NEEDLE

Code	Description	Drain size	Q.ty box	Drain type
MS/118V/Ch	Closed and continuous low vacuum drain set with collection bag 500 ml, connecting tube and silicone straight drain	From Ch 12 to Ch 30	30	
MS/119V/size	Closed and continuous low vacuum drain set with collection bag 500 ml, connecting tube and silicone flat drain 3/4 perforated			
MS/119VA/size	Same as above but with needle	7 mm x 20 cm 10 mm x 20 cm 13 mm x 20 cm	30	
MS/119VT/size	Closed and continuous low vacuum drain set with collection bag 500 ml, connecting tube and silicone flat drain 4/4 perforated			
MS/119VTA/size	Same as above but with needle			
MS/119V/size	Closed and continuous low vacuum drain set with collection bag 500 ml, connecting tube and silicone S-Line channelled drain, round or oval shaped	3, 5, 6, 7 mm (round type) 3.5 mm x 8 mm, 4 mm x 11 mm (flat type)	30	
MS/119VA/size	Same as above but with needle			



Note: replace "Ch" or "size" in the code with the drain size.

## REPLACEMENT BAG

Code	DescriptiON	Q.ty/box	
MS/115V/500	Spare collection bag 500 ml, with antibacterial air vent, non return valve	100	

## PVC DRAIN

- Various lengths and perforations
- ORX line
- With or without stainless steel needle
- Disposable, double packaged, sterile



Code	Description	Lenght	Size	Perforation	Q.ty/box
MS/106/50/Ch	PVC drain	50 cm	Ch. 06 - 18	7 cm	100
MS/107/50/Ch	PVC drain	50 cm	Ch. 06 - 18	15 cm	100
MS/107/80/Ch	PVC drain	80 cm	Ch. 06 - 18	15 cm	100
MS/108/50/Ch	PVC drain with needle	50 cm	Ch. 06 - 18	7 cm	50
MS/109/50/Ch	PVC drain with needle	50 cm	Ch. 06 - 18	15 cm	50
MS/109/80/Ch	PVC drain with needle	80 cm	Ch. 06 - 18	15 cm	50

Note: replace "Ch" in the code with the drain size (Ch. 06, 08, 10, 12, 14, 16 or 18).

