


Intraoperative suction tubes and cannulas



CANNULAS

- Various tips and handles available
- Complete range for intraoperative suction
- Q.ty per box: 100 pcs.

Standard cannula



Bubble cannula



Ergonomic cannula



Code	Description	Length	Vacuum control	Tip
MS/203	Bubble cannula with atraumatic soft tip	26 cm	no	
MS/203F	Bubble cannula with atraumatic soft tip	26 cm	yes	
MS/206	Ergonomic cannula (sump, central hole)	25 cm	no	
MS/206F	Ergonomic cannula (sump, central hole)	25 cm	yes	
MS/207	Bubble cannula with medium tip (Ch 18, six holes)	26 cm	no	
MS/207F	Bubble cannula with medium tip (Ch 18, six holes)	26 cm	yes	
MS/208	Bubble cannula with big tip (Ch 22, six holes)	26 cm	no	
MS/208F	Bubble cannula with big tip (Ch 22, six holes)	26 cm	yes	
MS/209	Bubble cannula with small tip (Ch 12, two holes)	26 cm	no	
MS/209F	Bubble cannula with small tip (Ch 12, two holes)	26 cm	yes	



CANNULAS WITH CONNECTING TUBE

- Set ready to use
- Tube available in various lengths
- Q.ty per box: 60 pcs. (for tube length 300 cm)

Code	Description	Length	Vacuum control	Tip
MS/203/300	Bubble cannula with atraumatic soft tip and connection tube 300 cm with female connector	26 cm	no	
MS/203F/300	Bubble cannula with atraumatic soft tip and connection tube 300 cm with female connector	26 cm	yes	
MS/206/300	Ergonomic cannula (sump, central hole) and connection tube 300 cm with female connector	25 cm	no	
MS/206F/300	Ergonomic cannula (sump, central hole) and connection tube 300 cm with female connector	25 cm	yes	
MS/207/300	Bubble cannula with medium tip (Ch. 18, six holes) and connection tube 300 cm with female connector	26 cm	no	
MS/207F/300	Bubble cannula with medium tip (Ch. 18, six holes) and connection tube 300 cm with female connector	26 cm	yes	
MS/208/300	Bubble cannula with big tip (Ch. 22, six holes) and connection tube 300 cm with female connector	26 cm	no	
MS/208F/300	Bubble cannula with big tip (Ch. 22, six holes) and connection tube 300 cm with female connector	26 cm	yes	
MS/209/300	Bubble cannula with small tip (Ch. 12, two holes) and connection tube 300 cm with female connector	26 cm	no	
MS/209F/300	Bubble cannula with small tip (Ch. 12, two holes) and connection tube 300 cm with female connector	26 cm	yes	

CONNECTION TUBES

Chimed manufactures a complete range of sterile surgical suction tubes, available in various length and diameters

Code	Description	Tube diameter	Connectors
MS/220/xxx	PVC connection tube length min 150 cm, max 600 cm	6,0 x 9,0 mm	
MS/230/xxx	Flex female/female connector	7,0 x10,0 mm	
MS/221/xxx	PVC connection tube length min 150 cm, max 600 cm	6,0 x 9,0 mm	
MS/231/xxx	Flex female/female connector with male adapter	7,0 x10,0 mm	

Tube length	Tube diameter	Q.ty/box
Up to 200 cm	6,0 x 9,0 mm or 7,0 x 10,0 mm	80
From 201 to 300 cm	6,0 x 9,0 mm	70
	7,0 x 10,0 mm	60
From 301 to 400 cm	6,0 x 9,0 mm	60
	7,0 x 10,0 mm	50
From 401 to 500 cm	6,0 x 9,0 mm	50
	7,0 x 10,0 mm	

