



Silicone drains and reservoirs



STRAIGHT DRAIN


- 100% silicone straight drain, 50 cm long
- With ORX line all along the tube
- 5 atraumatic side eyes at distal end and open tip
- Disposable, double packaged, sterile



Code	Description	Size	Length	Q.ty/box
MS/00013/Ch	Silicone straight drain	Ch 09, 12, 15, 16, 18, 20, 21, 22, 24, 27, 28, 30, 32, 33, 36	50 cm	10

FLAT DRAIN


- 100% silicone flat drain, with inner grooves
- Flat part completely ORX, $\frac{3}{4}$ or totally perforated
- Available with or without needle
- Disposable, double packaged, sterile



Code	Description	Size	Perforation	Q.ty/box
MS/00040/7x20	Silicone flat drain	3x7x200 mm	3/4	10
MS/00040/10x20	Silicone flat drain	4x10x200 mm	3/4	10
MS/00040/13x20	Silicone flat drain	5x13x200 mm	3/4	10
MS/4020/7x20	Silicone flat drain	3x7x200 mm	4/4	10
MS/4020/10x20	Silicone flat drain	4x10x200 mm	4/4	10
MS/4020/13x20	Silicone flat drain	5x13x200 mm	4/4	10
MS/00040A/7x20	Silicone flat drain with needle	3x7x200 mm	3/4	10
MS/00040A/10x20	Silicone flat drain with needle	4x10x200 mm	3/4	10
MS/00040A/13x20	Silicone flat drain with needle	5x13x200 mm	3/4	10
MS/4020A/7x20	Silicone flat drain with needle	3x7x200 mm	4/4	10
MS/4020A/10x20	Silicone flat drain with needle	4x10x200 mm	4/4	10
MS/4020A/13x20	Silicone flat drain with needle	5x13x200 mm	4/4	10

“T” DRAIN (Kehr drain)


- 100% silicone “T” drain
- Total length 70 cm.
- Cross line perforated or non perforated
- Disposable, double packaged, sterile



Code	Description	Size	Cross line	Q.ty/box
MS/00014/Ch	Silicone “T” drain	Ch. 10, 12, 14, 16, 18, 20, 22, 24	non perforated	10
MS/00015/Ch	Silicone “T” drain	Ch. 10, 12, 14, 16, 18, 20, 22 24	perforated	10

PENROSE DRAIN


- 100% silicone Penrose drain
- Length 30 cm o 50 cm
- Corrugated internal walls
- With ORX line all along the tube
- Disposable, double packaged, sterile



Code	Description	Size	Length	Q.ty/box
MS/00012/mm	Silicone Penrose drain	6, 8, 10, 12 mm	30 cm	25
MS/00012/50/mm	Silicone Penrose drain	6, 8, 10, 12 mm	50 cm	10

MULTILUMEN DRAIN (16 lumens)


- 100% silicone drain made up of 16 microlumens
- Length 30 or 40 cm
- Completely ORX
- Disposable, double packaged, sterile



Code	Description	Length	Q.ty/box
MS/00061	Silicone 16 multilumen drain	30 cm	10
MS/00061/40	Silicone 16 multilumen drain	40 cm	10

S-LINE CHANNELLED DRAIN


- 100% silicone drain with multiple narrow ducts
- Channelled part length 20 cm, ORX
- Tubular part length 60 cm, with or without needle
- Disposable, double packaged, sterile



Code	Description	Diameter	Design	Q.ty/box
MS/00084/03	Silicone S-Line drain	3 mm	round	10
MS/00084/05	Silicone S-Line drain	5 mm	round	10
MS/00084/06	Silicone S-Line drain	6 mm	round	10
MS/00084/07	Silicone S-Line drain	7 mm	round	10
MS/00084/3.5x8	Silicone S-Line drain	3,5 x 8 mm	flat	10
MS/00084/4x11	Silicone S-Line drain	4 x 11 mm	flat	10
MS/00085/03	Silicone S-Line drain with needle	3 mm	round	10
MS/00085/05	Silicone S-Line drain with needle	5 mm	round	10
MS/00085/06	Silicone S-Line drain with needle	6 mm	round	10
MS/00085/07	Silicone S-Line drain with needle	7 mm	round	10
MS/00085/3.5x8	Silicone S-Line drain with needle	3,5 x 8 mm	flat	10
MS/00085/4x11	Silicone S-Line drain with needle	4 x 11 mm	flat	10

RESERVOIRS

- 100% silicone oval shaped reservoir
- Capacities 100, 150, 200 or 400 ml
- Closing cap, non return valve
- Disposable, double packaged, sterile



Code	Description	Capacity	Drain connection	Q.ty/box
MS/00041/100	Silicone reservoir	100 ml	single	10
MS/00041/150	Silicone reservoir	150 ml	single	10
MS/00041/200	Silicone reservoir	200 ml	double	10
MS/00041/400	Silicone reservoir	400 ml	double	10

