


Passive drainage



Description

Each set includes:

- 1 collection bag, capacity 500 or 2000 ml, with non return valve, “T” cross outlet and bag hanger;
- 1 connection tube length 90 cm with connector for drain;
- 1 silicone straight drain, ORX,

Use

- Post-operative drainage with gravitational collection of liquids

Advantages

- Each device is supplied sterile double packaged.
- The closed system assures aseptic conditions avoiding external contamination.
- Double packaging for aseptic use in operating theatre.
- Transparent materials allow easy examination of the drained liquid.
- Spare parts available for cost reduction.

Code	Description	Drain type and sizes	Bag capacity	Q.ty/box
MS/120/Ch	Closed passive drain set	Straight, 50 cm long Ch. 12, 15, 18, 21, 24, 27, 30, 33, 36	2.000 ml	50
MS/121/Ch	Closed passive drain set	Straight, 50 cm long Ch. 12, 15, 18, 21, 24, 27, 30, 33, 36	500 ml	50

Note: replace “Ch” in the code with the drain size

